



AESCULAP® POWER SYSTEMS

CATALOGUE 2019

ELAN 4 and Acculan 4

CONTENT

ELAN 4

- ELAN 4 electro
- ELAN 4 air
- ELAN 4 Accessories

ACCULAN 4

- Acculan 4 Et Acculan 4 Small Drill
- Acculan 4 Electro
- Acculan 4 Accessories

PROCESSING

- Sterilit Power Systems
- Processing

OTHERS

- Tool Attachments Overview
 - Signs Et Symbols
 - Brochures
 - Index
-

AESCULAP® ELAN 4

PRECISION. PERFORMANCE. PLEASURE.

ELAN 4

ELAN 4 is the next generation of pneumatic and electric power systems for neuro and spine departments. With a new technical concept ELAN 4 addresses important requirements of the surgeon and all parties involved in handling power systems in the hospital. ELAN 4 electro and ELAN 4 air are two systems with one handling philosophy.

- Convenient and reliable grip
- Powerful, precise and reliable work
- Easy to set up and simple to use



DESIGN
AWARD
2016

The Winner of the
iF DESIGN AWARD 2016
Category Medicine/Healthcare

AESCULAP® ELAN 4

ELAN 4 electro

ELECTRIC HIGHSPEED SYSTEM



ELAN 4 electro

BASIC COMPONENTS

GA800

ELAN 4 electro CONTROL UNIT

with irrigation pump; touch screen;
power cord not included
pump flow 0–65 ml/min

100–240 V AC; 50/60 Hz;
330 x 330 x 200 mm; 9,500 g

consisting parts:
GD412804 bottle holder



GA808

ELAN 4 electro FOOT CONTROL

protection type IPX8; anaesthetic
proof class AP

5.0 m;
245 x 144 x 76 mm; 2,240 g



GA810

ELAN 4 electro WIRELESS FOOT CONTROL

for use with ELAN 4 electro applied
parts

271 x 266 x 85 mm; 3,800 g

Requires control unit software V3.0.



GA806

GB073R

ELAN 4 electro MOTOR CABLE

for use with ELAN 4 electro applied
parts

4.0 m; 270 g

consisting parts:
GB600860 Sterilit spray adaptor
ELAN 4 electro



AESCULAP® ELAN 4

ELAN 4 electro

BASIC COMPONENTS, PERFORATOR DRIVER & INTRA MOTOR

GA804

GB073R

ELAN 4 electro MOTOR CABLE WITH HAND SWITCH

for use with ELAN 4 electro applied
parts

4.0 m; 370 g

consisting parts:
GB600860 Sterilit spray adaptor
ELAN 4 electro



GA805

GB073R

ELAN 4 electro MOTOR CABLE WITH HAND CONTROL

for use with ELAN 4 electro applied
parts

4.0 m; 370 g

consisting parts:
GB600860 Sterilit spray adaptor
ELAN 4 electro



GA822

GB076R

ELAN 4 electro PERFORATOR DRIVER

for use with cranial perforators with
Hudson shank

180 W; 0 - 1,200 rpm; 250 Ncm;
166 x 35 mm; 410 g



GA824

GB077R

ELAN 4 electro LOWSPEED MOTOR

with intra coupling

180 W; 0 - 40,000* rpm;
10 Ncm; 123 x 23 mm; 150 g

* Requires control unit software V3.0.
With lower versions the motor runs
20,000 rpm.



ELAN 4 electro




CRANIOTOME & ATTACHMENTS

GA849

GB085R

ELAN 4 electro CRANIOTOME AND MULTIFUNCTION HANDPIECE

2-ring coding; for use with dura guards and holding sleeves GB941R – GB947R; spray nozzle GA259SU (not included)

140 W;  ; 0 – 80,000 rpm; 2.2 Ncm; 156 x 17 mm;  121 g



GB941R

GB719R

ELAN 4 DURA GUARD PEDIATRIC

fixed; 1 dot coding; for craniotome and multifunction handpieces GA749, GA849



GB942R

GB719R

ELAN 4 DURA GUARD STANDARD

fixed; 2 dot coding; for craniotome and multifunction handpieces GA749, GA849



GB943R

GB719R

ELAN 4 DURA GUARD LONG

fixed; 3 dot coding; for craniotome and multifunction handpieces GA749, GA849




GB945R

GB719R

ELAN 4 HOLDING SLEEVE

for craniotome and multifunction handpieces GA749, GA849



GB947R

GB719R

ELAN 4 DURA GUARD STANDARD

turnable; 2 dot coding; for craniotome and multifunction handpieces GA749, GA849



AESCULAP® ELAN 4




ELAN 4 electro

1-RING HIGHSPEED HANDPIECES

GA861

GB083R

ELAN 4 electro 1-RING HANDPIECE L4




1-ring coding;
spray nozzle GA261SU (not included)
140 W;  ; 0 – 80,000 rpm; 2.2 Ncm;
152 x 16 mm;  80 g



GA862

GB083R

ELAN 4 electro 1-RING HANDPIECE L7




1-ring coding;
spray nozzle GA262SU (not included)
140 W;  ; 0 – 80,000 rpm; 2.2 Ncm;
177 x 16 mm;  83 g



GA863

GB083R

ELAN 4 electro 1-RING HANDPIECE L10




1-ring coding;
spray nozzle GA263SU (not included)
140 W;  ; 0 – 80,000 rpm; 2.2 Ncm;
206 x 16 mm;  87 g



GA864

GB083R

ELAN 4 electro 1-RING HANDPIECE L13

1-ring coding;
spray nozzle GA264SU (not included)
140 W;  ; 0 – 80,000 rpm; 2.2 Ncm;
237 x 16 mm;  91 g






ELAN 4 electro

2-RING HIGHSPEED HANDPIECES

GA865

GB083R

ELAN 4 electro 2-RING HANDPIECE L7




2-ring coding;
spray nozzle GA265SU (not included)
140 W;  ; 0–80,000 rpm; 2.2 Ncm;
176 x 16 mm;  89 g




GA866

GB083R

ELAN 4 electro 2-RING HANDPIECE L10




2-ring coding;
spray nozzle GA266SU (not included)
140 W;  ; 0–80,000 rpm; 2.2 Ncm;
206 x 16 mm;  95 g




GA867

GB083R

ELAN 4 electro 2-RING HANDPIECE L13




2-ring coding;
spray nozzle GA267SU (not included)
140 W;  ; 0–80,000 rpm; 2.2 Ncm;
236 x 16 mm;  100 g



GA868

GB083R

ELAN 4 electro 2-RING HANDPIECE L17




2-ring coding;
spray nozzle GA268SU (not included)
140 W;  ; 0–80,000 rpm; 2.2 Ncm;
276 x 16 mm;  108 g



GA869

GB083R

ELAN 4 electro 2-RING HANDPIECE L22

2-ring coding;
spray nozzle GA269SU (not included)
140 W;  ; 0–80,000 rpm; 2.2 Ncm;
326 x 16 mm;  118 g



AESCULAP® ELAN 4

ELAN 4 electro

MIS HANDPIECE & HANDPIECE SHAFTS



GA860**GB084R****ELAN 4 electro MIS HANDPIECE**for use with handpiece shafts
GB920R - GB926R140 W; ; 0 - 80,000 rpm;
2.2 Ncm; 118 x 17 mm;  100 g**GB920R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L10 STRAIGHT**

spray nozzle GB796SU (not included)

**GB921R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L10 CURVED**

spray nozzle GB796SU (not included)

**GB922R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L10 HEAVY CURVED**

spray nozzle GB796SU (not included)

**GB925R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L13 STRAIGHT**

spray nozzle GB797SU (not included)

**GB926R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L13 CURVED**

spray nozzle GB797SU (not included)

**GB927R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L13 HEAVY CURVED**

spray nozzle GB797SU (not included)



AESCULAP® ELAN 4


ELAN 4 electro

SAWS



GA836**GB086R****ELAN 4 electro MICRO SAGITTAL SAW**

spray nozzle GA258SU (not included)

180 W; 0 - 20,000 rpm; 161 x 23 mm;
 \angle 7°; 190 g**GA837****GB087R****ELAN 4 electro MICRO RECIPROCATING SAW**

spray nozzle GA258SU (not included)

180 W; 0 - 20,000 rpm; 184 x 23 mm;
 \updownarrow 1.3 mm; 198 g**GA831****GB088R****ELAN 4 electro SAGITTAL SAW**


spray nozzle GA258SU (not included)

180 W; 0 - 20,000 rpm; 129 x 23 mm;
 \angle 4.5°; 220 g**GA832****GB087R****ELAN 4 electro RECIPROCATING SAW**

spray nozzle GA258SU (not included)

180 W; 0 - 20,000 rpm; 194 x 23 mm;
 \updownarrow 2.5 mm; 238 g**GA833****GB088R****ELAN 4 electro TRANSVERSAL SAW**

spray nozzle GA258SU (not included)

180 W; $\swarrow \searrow$; 0 - 20,000 rpm; 10 Ncm;
 130 x 23 mm; \angle 10°; 185 g

AESCULAP® ELAN 4

ELAN 4 electro

DRILL & ATTACHMENTS

GA844

GB482R

ELAN 4 electro DRILL

with oscillating and tapping mode

200 W;  3.3 mm;

130 x 27 x 122 mm;  800 g

consisting parts:

GA344244 K-wire protection guide

GA344211 rinsing adaptor

GB600880 Sterilit spray adaptor

Requires control unit software V3.0.





GB870R

GB716R

INTRAMEDULLARY REAMER ATTACHMENT

with Jacobs chuck 0.5 – 7.4 mm,
not for use with acetabulum reamers

0 – 320 rpm; 6.5 Nm;  3.2 mm;

105 x 29 mm;  230 g

consisting parts:

key GA031R





GB871R

GB716R

INTRAMEDULLARY REAMER ATTACHMENT

with large AO chuck,
not for use with acetabulum reamers

0 – 320 rpm; 6.5 Nm;  3.2 mm;

88 x 29 mm;  170 g





GB872R

GB716R

INTRAMEDULLARY REAMER ATTACHMENT

with Hudson/Zimmer chuck,
not for use with acetabulum reamers

0 – 320 rpm; 6.5 Nm;  3.2 mm;

86 x 25 mm;  150 g



ELAN 4 electro ATTACHMENTS

GB881R

GB716R

DRILL ATTACHMENT

with Jacobs chuck 0.5 – 7.4 mm

0 – 1,250 rpm; 2 Nm;  3.2 mm;
89 x 29 mm;  163 g

consisting parts:
key GA031R




GB882R

GB716R

DRILL ATTACHMENT

with Jacobs chuck 0.5 – 4 mm

0 – 1,250 rpm; 2 Nm;  3.2 mm;
76 x 22 mm;  100 g

consisting parts:
key GA062R





GB884R

GB716R

DRILL ATTACHMENT

with small A0 chuck

0 – 1,250 rpm; 2 Nm;  2.8 mm;
60 x 22 mm;  66 g





GB886R

GB716R

DRILL ATTACHMENT

with AESCULAP® hexagonal chuck

0 – 1,250 rpm; 2 Nm;  3.2 mm;
57 x 24 mm;  67 g





GB887R

GB716R

DRILL ATTACHMENT

with Dental shank

0 – 1,250 rpm; 2 Nm;  1.8 mm;
56 x 22 mm;  50 g



AESCULAP® ELAN 4

ELAN 4 electro

ATTACHMENTS

GB878R

GB716R

ADAPTOR

for Synthes radiolucent drive 511.300

0 – 1,250 rpm; 2 Nm;
111 x 27 mm; 181 g



GB894R

GB716R

K-WIRE ATTACHMENT

for K-wires from 0.6 – 3.2 mm

0 – 1,250 rpm; 3.2 mm;
92 x 27 x 124 mm; 179 g



GB891R

GB716R

SAGITTAL SAW ATTACHMENT

for "L-sagittal" sawblades with a
maximum usable length of 50 mm

0 – 15,750 rpm; 74 x 33.5 mm;
4.5°; 179 g



GB892R

GB716R

RECIPROCATING SAW ATTACHMENT



for reciprocating sawblades

0 – 15,750 rpm; 135 x 24.5 mm;
2.5 mm; 166 g





GB896R**GB716R****SCREW ATTACHMENT**

with small AO chuck

0 - 165 rpm; 12 Nm;  1.7 mm;
95 x 28 x 129.5 mm;  235 g**GB897R****GB716R****SCREW ATTACHMENT**

with AESCULAP® hexagonal chuck

0 - 165 rpm; 12 Nm;  1.7 mm;
95 x 28 x 129.5 mm;  231 g

AESCULAP® ELAN 4

ELAN 4 air

PNEUMATIC HIGHSPEED SYSTEM



ELAN 4 air WALL HOSE & ADAPTOR

GA702R

ELAN 4 air WALL HOSE

for use with ELAN 4 air wall adaptor
GA710R - GA715R

 3.0 m

GA703R

ELAN 4 air WALL HOSE

for use with ELAN 4 air wall adaptor
GA710R - GA715R

 5.0 m



GA710R

ELAN 4 air WALL ADAPTOR

AESCULAP® Draeger

for use with ELAN 4 air hoses
GA702R, GA703R, GA706R, GA707R



GA712R

ELAN 4 air WALL ADAPTOR

Schrader with diffuser

for use with ELAN 4 air hoses
GA702R, GA703R, GA706R, GA707R



GA713R

ELAN 4 air WALL ADAPTOR

DIN with diffuser

for use with ELAN 4 air hoses
GA702R, GA703R, GA706R, GA707R



GA715R

ELAN 4 air WALL ADAPTOR

Synthes

for use with ELAN 4 air hoses
GA702R, GA703R, GA706R, GA707R



AESCULAP® ELAN 4

ELAN 4 air

FOOT CONTROL, MOTOR HOSE & PERFORATOR DRIVER

GA708

ELAN 4 air FOOT CONTROL

for use with ELAN 4 air wall adaptor GA710R – GA715R; ELAN 4 air wall hose GA702R, GA703R; ELAN 4 air motor hose GA705R – GA707R



GA705R

GB053R

ELAN 4 air MOTOR HOSE

for use with ELAN 4 air applied parts and ELAN 4 air foot control GA708

3.0 m

consisting parts:
GB600850 Sterilit spray adaptor
ELAN 4 air



GA706R

GB053R

ELAN 4 air MOTOR HOSE

with hand control; for use with ELAN 4 air applied parts and ELAN 4 wall adaptor GA710R – GA715R

3.0 m

consisting parts:
GB600850 Sterilit spray adaptor
ELAN 4 air



GA707R


GB053R

ELAN 4 air MOTOR HOSE

with hand control; for use with ELAN 4 air applied parts and ELAN 4 wall adaptor GA710R – GA715R

5.0 m

consisting parts:
GB600850 Sterilit spray adaptor
ELAN 4 air



GA722

GB056R

ELAN 4 air PERFORATOR DRIVER

for use with cranial perforators with Hudson shank

100 W; 0 – 1,200 rpm; 400 Ncm;
155 x 35 mm; 427 g



ELAN 4 air

CRANIOTOME & ATTACHMENTS

GA749

GB065R

ELAN 4 air CRANIOTOME AND MULTIFUNCTION HANDPIECE

2-ring coding; for use with dura guards and holding sleeves GB941R – GB947R

80 W; 0 – 80,000 rpm; 3.5 Ncm;
128 x 18 mm; 108 g



GB941R

GB719R

ELAN 4 DURA GUARD PEDIATRIC

fixed; 1 dot coding;
for craniotome and multifunction
handpieces GA749, GA849



GB942R

GB719R

ELAN 4 DURA GUARD STANDARD

fixed; 2 dot coding;
for craniotome and multifunction
handpieces GA749, GA849



GB943R

GB719R

ELAN 4 DURA GUARD LONG

fixed; 3 dot coding;
for craniotome and multifunction
handpieces GA749, GA849




GB945R

GB719R

ELAN 4 HOLDING SLEEVE

for craniotome and multifunction
handpieces GA749, GA849



GB947R

GB719R

ELAN 4 DURA GUARD STANDARD

turnable; 2 dot coding;
for craniotome and multifunction
handpieces GA749, GA849



AESCULAP® ELAN 4

ELAN 4 air



1-RING HIGHSPEED HANDPIECES

GA761

GB063R

ELAN 4 air 1-RING HANDPIECE L4

1-ring coding

80 W;  0-80,000 rpm; 3.5 Ncm;
125 x 18 mm;  68 g






GA762

GB063R

ELAN 4 air 1-RING HANDPIECE L7

1-ring coding

80 W;  0-80,000 rpm; 3.5 Ncm;
150 x 18 mm;  71 g






GA763

GB063R

ELAN 4 air 1-RING HANDPIECE L10

1-ring coding

80 W;  0-80,000 rpm; 3.5 Ncm;
180 x 18 mm;  75 g




GA764

GB063R

ELAN 4 air 1-RING HANDPIECE L13

1-ring coding

80 W;  0-80,000 rpm; 3.5 Ncm;
210 x 18 mm;  79 g



ELAN 4 air

2-RING HIGHSPEED HANDPIECES


GA765

GB063R

ELAN 4 air 2-RING HANDPIECE L7

2-ring coding

80 W; 0-80,000 rpm; 3.5 Ncm;
150 x 18 mm; 76 g




GA766

GB063R

ELAN 4 air 2-RING HANDPIECE L10

2-ring coding

80 W; 0-80,000 rpm; 3.5 Ncm;
180 x 18 mm; 82 g



GA767

GB063R

ELAN 4 air 2-RING HANDPIECE L13

2-ring coding

80 W; 0-80,000 rpm; 3.5 Ncm;
210 x 18 mm; 88 g



GA768

GB063R

ELAN 4 air 2-RING HANDPIECE L17

2-ring coding

80 W; 0-80,000 rpm; 3.5 Ncm;
250 x 18 mm; 95 g



GA769

GB063R

ELAN 4 air 2-RING HANDPIECE L22

2-ring coding

80 W; 0-80,000 rpm; 3.5 Ncm;
300 x 18 mm; 105 g



AESCULAP® ELAN 4

ELAN 4 air

MIS HANDPIECE & HANDPIECE SHAFTS



GA760**GB064R****ELAN 4 air MIS HANDPIECE**for use with handpiece shafts
GB920R - GB927R80 W; 0 - 80,000 rpm;
3.5 Ncm; 89 x 18 mm; 85 g**GB920R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L10 STRAIGHT****GB921R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L10 CURVED****GB922R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L10 HEAVY CURVED****GB925R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L13 STRAIGHT****GB926R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L13 CURVED****GB927R****GB723R****ELAN 4 MIS HANDPIECE SHAFT
L13 HEAVY CURVED**

AESCULAP® ELAN 4

ELAN 4 air

SAWS & INTRA MOTOR



GA736**GB066R****ELAN 4 air MICRO SAGITTAL SAW**

60 W; 0 – 18,000 rpm; 150 x 23 mm;
 \angle 7°; \updownarrow 214 g

**GA737****GB067R****ELAN 4 air MICRO RECIPROCATING SAW**

75 W; 0 – 21,000 rpm; 167 x 23 mm;
 \updownarrow 1.3 mm; \updownarrow 215 g

**GA731****GB068R****ELAN 4 air SAGITTAL SAW**

75 W; 0 – 21,000 rpm; 116 x 23 mm;
 \angle 4.5°; \updownarrow 240 g

**GA732****GB067R****ELAN 4 air RECIPROCATING SAW**

75 W; 0 – 21,000 rpm; 118 x 23 mm;
 \updownarrow 2.5 mm; \updownarrow 250 g

**GA733****GB068R****ELAN 4 air TRANSVERSAL SAW**

75 W; \searrow ; 0 – 21,000 rpm; 7.6 Ncm;
 115 x 23 mm; \angle 10°; \updownarrow 200 g

**GA724****GB057R****ELAN 4 air LOWSPEED MOTOR**

with intra coupling

80 W; \searrow ; 0 – 21,000 rpm; 7.6 Ncm;
 93 x 23 mm; \updownarrow 130 g



AESCLAP[®] ELAN 4

ELAN 4 ACCESSORIES

FOR ELAN 4 air & ELAN 4 electro

GB861P

GB727R

ELAN 4 LENGTH REDUCTION SLEEVE

for use with 1-ring handpieces



GB862P

GB727R

ELAN 4 MIS SLEEVE

for use with 1-ring handpieces



GB863P

GB727R

ELAN 4 LENGTH REDUCTION SLEEVE

for use with 2-ring handpieces



GB866R

ELAN 4 LATERAL BURR GUARD

for use with 1-ring handpieces and
burrs up to 4.5 mm diameter



TE780

POWER CORD

Europe, 1.5 m

TE730

POWER CORD

Europe, 5.0 m

TE734

POWER CORD

UK, 5.0 m

TE735


POWER CORD

USA, Canada, Japan, 3.5 m



GA395SUsterile  

ELAN 4 electro DISPOSABLE TUBE SET

**ELAN 4 SPRAY NOZZLES**sterile **GA259SU**

for craniotome attachment

GA261SU

for 1-ring handpiece L4

GA262SU

for 1-ring handpiece L7

GA263SU

for 1-ring handpiece L10

GA264SU

for 1-ring handpiece L13

GA265SU

for 2-ring handpiece L7

GA266SU

for 2-ring handpiece L10

GA267SU

for 2-ring handpiece L13

GA268SU

for 2-ring handpiece L17

GA269SU

for 2-ring handpiece L22

GB796SU

for MIS handpiece shafts L10

GB797SU

for MIS handpiece shafts L13

GA258SU

for saws



AESCULAP® Acculan 4

LOOK. FEEL. EXPERIENCE.

Acculan 4 the battery-operated power system for orthopaedics, traumatology, cardiothoracic surgery and dermatology

- Powerful
- Longevity – Exclusion of battery and electronic control from sterilization process
- Easy handling – Titanium housing allows mechanical alkaline reprocessing

ACCULAN 4



AESCULAP® Acculan 4

DRILL & ATTACHMENTS

GA330

GB495R

Acculan 4 DRILL AND REAMER

250 W;   4.0 mm;
207 x 121 x 55 mm;  1,250 g

consisting parts:

GA678 sterile funnel

GA672207 K-wire protection sleeve

GA672211 rinsing adaptor

GB600840 Sterilit spray adaptor

Acculan

Lid GA675 and NiMH battery GA676 to be ordered separately.

Measurements and weight are including battery and lid.





GB620R

GB496R

Acculan 4 DRILL ATTACHMENT

with Jacobs chuck keyless 0.5 – 7.4 mm

0 – 1,000 rpm; 5.0 Nm;  4.0 mm;
133 x 35 mm;  397 g




GB621R

GB496R

Acculan 4 DRILL ATTACHMENT

with Jacobs chuck 0.5 – 7.4 mm

0 – 1,000 rpm; 5.0 Nm;  4.0 mm;
118 x 35 mm;  300 g

consisting parts:

key GA031R





GB623R

GB496R

Acculan 4 DRILL ATTACHMENT

with small AO chuck

0 – 1,250 rpm; 4.0 Nm;  2.8 mm;
90 x 35 mm;  210 g





ATTACHMENTS

GB639R

GB496R

Acculan 4 DRILL ATTACHMENT

with Jacobs chuck 0.5 – 4.0 mm

0 – 1,250 rpm; 4.0 Nm;  4.0 mm;
105 x 35 mm;  243 g

consisting parts:
key GA062R



GB645R

GB496R

Acculan DRILL ATTACHMENT

for Synthes radiolucent drive 511.300

0 – 900 rpm; 110 x 35 mm;  225 g





GB663R

GB496R

Acculan 4 DRILL ATTACHMENT

with trinkle AO chuck

0 – 1,000 rpm; 5.0 Nm;  4.0 mm;
93 x 35 mm;  217 g





GB664R

GB496R

Acculan 4 DRILL ATTACHMENT

with AESCULAP® hexagonal chuck

0 – 1,000 rpm; 5.0 Nm;  4.0 mm;
86 x 35 mm;  225 g





GB665R

GB496R

Acculan 4 DRILL ATTACHMENT

with Hudson/Zimmer chuck

0 – 1,000 rpm; 5.0 Nm;  4.0 mm;
100 x 35 mm;  230 g



AESCULAP® Acculan 4



ATTACHMENTS

GB657R

GB496R

Acculan 4 INTRAMEDULLAR ATTACHMENT

with large AO chuck

0–370 rpm; 12.0 Nm;  4.0 mm;
108 x 35 mm;  300 g





GB667R

GB496R

Acculan 4 REAMER ATTACHMENT

with Jacobs chuck 0.5–7.4 mm chuck

0–250 rpm; 19.0 Nm;  4.0 mm;
133 x 35 mm;  280 g

consisting parts:
key GA031R





GB668R

GB496R

Acculan 4 REAMER ATTACHMENT

with large AO chuck

0–250 rpm; 19.0 Nm;  4.0 mm;
108 x 35 mm;  300 g





GB669R

GB496R

Acculan 4 REAMER ATTACHMENT

with Hudson/Zimmer chuck

0–250 rpm; 19.0 Nm;  4.0 mm;
118 x 35 mm;  300 g





GB670R

GB496R

Acculan 4 REAMER ATTACHMENT

with Harris chuck

0–250 rpm; 19.0 Nm;  4.0 mm;
106 x 35 mm;  290 g



GB660R**GB496R****Acculan SAGITTAL SAW ATTACHMENT**

for Acculan rapid action saw blades
with a maximum useable length of
50 mm

«» 0 - 17,000 rpm; \angle 4°47';
89 x 36 x 44 mm; ■ 270 g

**GB641R****GB496R****Acculan K-WIRE ATTACHMENT**

for K-wires from 0.6 - 4.0 mm

0 - 1,250 rpm; 93 x 36 x 137 mm;
 ■ 285 g



AESCULAP® Acculan 4

OSCILLATING SAWS

GA331

GB495R

Acculan 4 OSCILLATING SAW

250 W; \leftrightarrow 0 – 13,000 rpm; \angle 4.1°;
211 x 177 x 55 mm; 1,600 g

consisting parts:
GA678 sterile funnel

Lid GA675 and NiMH battery GA676 to
be ordered separately.

Measurements and weight are including
battery and lid.



GA336

GB498R

Acculan 4 OSCILLATING SAW STRAIGHT

250 W; \leftrightarrow 0 – 13,000 rpm; \angle 4.1°;
290 x 55 x 84 mm; 1,480 g

consisting parts:
GA678 sterile funnel

Lid GA675 and NiMH battery GA676 to
be ordered separately.

Measurements and weight are including
battery and lid.



GA337

GB498R

Acculan 4 OSCILLATING STERNUM SAW

250 W; \leftrightarrow 0 – 13,000 rpm; \angle 4.1°;
290 x 55 x 84 mm; 1,480 g

consisting parts:
GA678 sterile funnel

Lid GA675 and NiMH battery GA676 to
be ordered separately.

Measurements and weight are including
battery and lid.



RECIPROCATING SAW

GA334

GB495R

Acculan 4 RECIPROCATING SAW

250 W; ⇄ 0 – 15,000 rpm; ↓ 3.2 mm;
207 x 180 x 55 mm; ⚖ 1,540 g

consisting parts:

GA678 sterile funnel

GB440215 cleaning plate

GA674890 cleaning adaptor

Lid GA675 and NiMH battery GA676 to
be ordered separately.

Measurements and weight are including
battery and lid.



GB436R

GB686R

PROTECTION GUIDE

for use with GA334 and saw blade
GC761R



AESCULAP® Acculan 4

DERMATOME

GA340

GB498R

Acculan 4 DERMATOME

width keyless adjustable;
thickness adjustable (1 / 10 mm steps);
cutting thickness: 0.2 - 1.2 mm;
cutting width: max. 78 mm;
suitable blades: GB228R
(to be ordered separately)

250 W; ⇄ 0 - 6,500 rpm; ⬆ 3.1 mm;
277 x 113 x 83 mm; ⚖ 1,650 g

consisting parts:
GA678 sterile funnel

Lid GA675 and small Li-Ion battery
GA346 to be ordered separately.

Measurements and weight are including
battery and lid.



GA341

GB498R

Acculan 4 DERMATOME 0.1 MM

width keyless adjustable;
thickness adjustable (1 / 10 mm steps);
cutting thickness: 0.1 - 1.1 mm;
cutting width: max. 78 mm;
suitable blades: GB228R
(to be ordered separately)

250 W; ⇄ 0 - 6,500 rpm; ⬆ 3.1 mm;
277 x 113 x 83 mm; ⚖ 1,650 g

consisting parts:
GA678 sterile funnel

Lid GA675 and small Li-Ion battery
GA346 to be ordered separately.

Measurements and weight are including
battery and lid.



BA720R

GB688R

SKIN GRAFT MESHER

for use in combination with skin graft
carriers BA721 (1:1.5) / BA722 (1:3) /
BA723 (1:6)

150 x 195 x 105 mm; ⚖ 4,600 g

consisting parts:
BA725R cutting roller
BA726R ratchet wrench



AESCULAP® Acculan 4

LOOK. FEEL. EXPERIENCE.

Acculan 4 small drill – the smallest member of the Acculan 4 family can be used for small fragment Osteosynthesis as well as for reamed intramedullary nailing.

ACCULAN 4



AESCULAP® Acculan 4

SMALL DRILL

DRILL & ATTACHMENTS

GA344

GB482R

Acculan 4 SMALL DRILL

with oscillating and tapping mode;
not for use with acetabulum reamers

85 W;  3.2 mm;
177 x 124 x 55 mm;  840 g

consisting parts:

GA678 sterile funnel

GA344244 K-wire protection guide

GA344211 rinsing adaptor

GB600880 Sterilit spray adaptor

Lid GA675 and small Li-Ion battery
GA346 to be ordered separately.

Measurements and weight are including
battery and lid.





GB870R

GB716R

INTRAMEDULLARY REAMER ATTACHMENT

with Jacobs chuck 0.5 - 7.4 mm,
not for use with acetabulum reamers

0 - 320 rpm; 5 Nm;  3.2 mm;
105 x 29 mm;  230 g

consisting parts:

key GA031R



GB871R

GB716R

INTRAMEDULLARY REAMER ATTACHMENT

with large AO chuck,
not for use with acetabulum reamers

0 - 320 rpm; 5 Nm;  3.2 mm;
88 x 29 mm;  170 g





GB872R

GB716R

INTRAMEDULLARY REAMER ATTACHMENT

with Hudson/Zimmer chuck,
not for use with acetabulum reamers

0 - 320 rpm; 5 Nm;  3.2 mm;
86 x 25 mm;  150 g





SMALL DRILL ATTACHMENTS

GB881R

GB716R

DRILL ATTACHMENT

with Jacobs chuck 0.5 – 7.4 mm

0 – 1,250 rpm; 1.5 Nm;  3.2 mm;
89 x 29 mm;  163 g

consisting parts:
key GA031R





GB882R

GB716R

DRILL ATTACHMENT

with Jacobs chuck 0.5 – 4 mm

0 – 1,250 rpm; 1.5 Nm;  3.2 mm;
76 x 22 mm;  100 g

consisting parts:
key GA062R





GB884R

GB716R

DRILL ATTACHMENT

with small A0 chuck

0 – 1,250 rpm; 1.5 Nm;  2.8 mm;
60 x 22 mm;  66 g





GB886R

GB716R

DRILL ATTACHMENT

with AESCULAP® hexagonal chuck

0 – 1,250 rpm; 1.5 Nm;  3.2 mm;
57 x 24 mm;  67 g





GB887R

GB716R

DRILL ATTACHMENT

with Dental shank

0 – 1,250 rpm; 1.5 Nm;  1.8 mm;
56 x 22 mm;  50 g



AESCULAP® Acculan 4

SMALL DRILL ATTACHMENTS

GB878R

GB716R

ADAPTOR

for Synthes radiolucent drive 511.300

0 – 1,250 rpm; 1.5 Nm;
111 x 27 mm; 181 g



GB894R

GB716R

K-WIRE ATTACHMENT

for K-wires from 0.6 – 3.2 mm

0 – 1,250 rpm; 3.2 mm;
92 x 27 x 124 mm; 179 g



GB891R

GB716R

SAGITTAL SAW ATTACHMENT

for "L-sagittal" sawblades with a
maximum usable length of 50 mm

0 – 15,750 rpm; 74 x 33.5 mm;
4.5°; 179 g



GB892R

GB716R

RECIPROCATING SAW ATTACHMENT



for reciprocating saw blades

0 – 15,750 rpm; 135 x 24.5 mm;
166 g; 2.5 mm





GB896R**GB716R****SCREW ATTACHMENT**

with small AO chuck

0 - 165 rpm; 10 Nm;  1.7 mm;
95 x 28 x 129.5 mm;  235 g**GB897R****GB716R****SCREW ATTACHMENT**

with AESCULAP® hexagonal chuck

0 - 165 rpm; 10 Nm;  1.7 mm;
95 x 28 x 129.5 mm;  231 g

AESCULAP® Acculan 4

CHARGER & BATTERIES

GA320

Acculan 4 CHARGER

for Acculan Li-Ion and NiMH batteries GA346, GA676; 4 charging bays; automatic battery diagnosis function; power cord not included

100 – 240 V AC; 50 / 60 Hz;
410 x 120 x 160 mm; 2,300 g



GA346

Acculan 4 LI-ION BATTERY SHORT

only for GA340, GA341 and GA344;
with integrated electronic motor control 9.9 V; 10.89 Wh

120 x 43 x 50 mm; 265 g



GA676

Acculan NIMH BATTERY

only for GA330, GA331, GA334, GA336, GA337; with integrated electronic motor control 9.6 V; 1.95 Ah

140 x 43 x 50 mm; 425 g



GA679

GB489R

Acculan BATTERY REMOVAL DEVICE



GA675

GB487R

Acculan LID

90 g



GA678

GB488R

Acculan STERILE FUNNEL



AESCULAP® Acculan 4 Electro

LOOK. FEEL. EXPERIENCE.




CONSOLE & BATTERIES

GA877

Acculan Electro CONSOLE

power cord not included


100 – 240 V AC; 50 / 60 Hz;
255 x 217 x 90 mm; 2,930 g



GA878 GB486R

Acculan Electro CONNECTION CABLE

for use with Acculan 4 handpieces

 3.0 m; 410 g



GA875 GB485R

Acculan Electro LID



GA876

Acculan Electro ADAPTOR

120 x 49 x 42 mm; 165 g





GA678 GB488R

Acculan STERILE FUNNEL



AESCULAP® Acculan 4 ACCESSORIES

ADAPTORS, KEYS & POWER CORD

GB176R GB467R ZIMMER ADAPTOR tool: Zimmer motor: Hudson/Zimmer only in combination with GB630R, GB656R, GB665R, GB669R or GB872R  3.8 mm	GB628R GB690R SYNTHES DHS/DCS ADAPTOR tool: Synthes DHS/DCS motor: AO large  4.2 mm	GB184R GB690R HARRIS ADAPTOR tool: Harris motor: AO large  3.5 mm
XF457R GB690R DIN ADAPTOR tool: DIN motor: AO large  3.8 mm	GA031R KEY for large Jacobs chuck suitable for GB620R, GB621R, GB667R, GB870R, GB881R	GA062R GB690R KEY for small 4 mm Jacobs chuck suitable for GB639R, GB882R
TE780 POWER CORD Europe,  1.5 m		
TE730 POWER CORD Europe,  5.0 m		
TE734 POWER CORD UK,  5.0 m		
TE735 POWER CORD USA, Canada, Japan,  3.5 m		

AESCULAP® PROCESSING

Sterilit POWER SYSTEMS



AESCULAP® PROCESSING

Sterilit POWER SYSTEMS

GB600



Sterilit POWER SYSTEMS
oil spray; 300 ml



GB600850

Sterilit POWER SYSTEMS SPRAY
ADAPTOR
for ELAN 4 air



GB600860

Sterilit POWER SYSTEMS SPRAY
ADAPTOR
for ELAN 4 electro



GB600870

Sterilit POWER SYSTEMS SPRAY
ADAPTOR
for ELAN 4 MIS handpiece shafts



GB600840

Sterilit POWER SYSTEMS SPRAY
ADAPTOR
for Acculan 4



GB600880

Sterilit POWER SYSTEMS SPRAY
ADAPTOR
for Acculan 4 small drill and ELAN 4
electro drill



AESFULAP® MECHANICAL PROCESSING

- The solution for storage and mechanical reprocessing of our AESFULAP® Power Systems.
- For detailed information please refer to the Instruction for Use.



AESCU LAP[®] MECHANICAL PROCESSING

ECCOS SETS

ELAN 4

GB261R

ELAN 4 electro ECCOS SET

consisting parts:

GB073R

GB076R


GB083R

GB085R

GB719R

GB072R

JF212R



GB260R

ELAN 4 air ECCOS SET

consisting parts:

GB053R

GB056R


GB063R

GB065R

GB719R

GB072R

JF212R



GB262R

ELAN 4 electro ECCOS SET FOR DRILL GA844

consisting parts:

GB072R

GB482R


GB715R

GB716R

JG340




























JG342

JG344



ECCOS









ELAN 4 electro

GB073R ECCOS HOLDER for motor cable GA804, GA805, GA806 (set consisting of 2 pieces)  56 mm	GB076R ECCOS HOLDER for perforator driver GA822  56 mm	GB077R ECCOS HOLDER for intra motor GA824  56 mm
 	 	 
GB082R ECCOS HOLDER for 1 handpiece  56 mm	GB083R ECCOS HOLDER for 3 handpieces  56 mm	GB085R ECCOS HOLDER for craniotome GA849  56 mm
 	 	 
GB086R ECCOS HOLDER for micro sagittal saw GA836  56 mm	GB087R ECCOS HOLDER for reciprocating saws GA832, GA837  56 mm	GB088R ECCOS HOLDER for saws GA831, GA833  56 mm
 	 	 

AESULAP[®] MECHANICAL PROCESSING











ECCOS

ELAN 4 electro

GB482R	GB715R	GB716R
<p>ECCOS HOLDER</p> <p>for ELAN 4 electro drill GA844</p> <p> 106 mm</p>	<p>ECCOS HOLDER</p> <p>for 1 attachment GB870R, GB871R, GB872R, GB878R, GB881R, GB884R, GB892R, GB896R, GB897R</p> <p> 106 mm</p>	<p>ECCOS HOLDER</p> <p>for 3 attachments GB870R, GB871R, GB872R, GB878R, GB881R, GB884R, GB892R, GB896R, GB897R</p> <p> 106 mm</p>
	 	 

ECCOS & RINSING DEVICES

ELAN 4 electro

<div>GB084R</div> <div>ECCOS HOLDER</div> <div>for MIS handpiece GA860</div> <div>   </div>	<div>GB723R</div> <div>ECCOS HOLDER</div> <div>for MIS handpiece shafts GB920R - GB927R</div> <div>   </div>
<div>GB692R</div> <div>RINSING DEVICE</div> <div>for applied parts, incl. ceramic filter GB693R (not suitable for autoclaving)</div> <div>   </div>	<div>GB679R</div> <div>RINSING DEVICE</div> <div>for 4 MIS shafts, incl. ceramic filter GB693R (not suitable for autoclaving)</div> <div>   </div>
<div>GB698R</div> <div>RINSING ADAPTOR</div> <div>for applied parts (for manual cleaning)</div> <div>  </div>	<div>GB699R</div> <div>RINSING ADAPTOR</div> <div>for MIS shafts (for manual cleaning)</div> <div>  </div>

AESCLAP[®] MECHANICAL PROCESSING

ECCOS

ELAN 4 air

GB053R

ECCOS HOLDER

for motor hose GA705R, GA706R,
GA707R

 56 mm



GB056R

ECCOS HOLDER

for perforator driver GA722

 56 mm




GB057R

ECCOS HOLDER

for intra motor GA724

 56 mm




GB062R

ECCOS HOLDER

for 1 handpiece

 56 mm




GB063R

ECCOS HOLDER

for 3 handpieces

 56 mm



GB065R

ECCOS HOLDER

for craniotome GA749

 56 mm



GB066R

ECCOS HOLDER

for micro sagittal saw GA736

 56 mm



GB067R

ECCOS HOLDER

for reciprocating saws GA732, GA737

 56 mm



GB068R

ECCOS HOLDER

for saws GA731, GA733

 56 mm




ECCOS & RINSING DEVICES

ELAN 4 air

GB064R

ECCOS HOLDER

for MIS handpiece GA760



GB723R

ECCOS HOLDER


for 2 MIS shafts GB920R - GB927R



GB691R

RINSING DEVICE


for applied parts, incl. ceramic filter GB693R
(not suitable for autoclaving)



GB679R

RINSING DEVICE

for 4 MIS shafts, incl. ceramic filter GB693R
(not suitable for autoclaving)



GB697R

RINSING ADAPTOR

for applied parts
(for manual cleaning)



GB699R

RINSING ADAPTOR
















for MIS shafts
(for manual cleaning)



AESCU LAP[®] MECHANICAL PROCESSING

ECCOS

ELAN 4 ACCESSORIES

GB718R ECCOS HOLDER for 12 1- and 2-ring burrs	GB719R ECCOS HOLDER for 3 craniotome attachments	GB720R ECCOS UNIVERSAL HOLDER for saw blades
 	 	 
GB729R ECCOS HOLDER for 2 reciprocating saw shafts	GB730R ECCOS HOLDER for 2 transversal saw shafts	GB072R ECCOS CABLE/HOSE POSITIONING GUIDE for cable and hoses (set consisting of 4 pieces)
 	 	 
GB727R ECCOS HOLDER for handpiece holding-sleeves GB861P, GB862P, GB863P	GB572R ECCOS FIXATION for GB300R - GB308R	
 	 	

ECCOS SETS

Acculan 4

GB257R

ECCOS SET

for Acculan 4 small drill GA344

 106 mm

consisting parts:

GB482R

GB487R


GB715R

GB716R

JG340

JG342

JG344



GB243R

ECCOS SET

for 2 Acculan 4 machines

 106 mm

consisting parts:

GB481R

GB487R


GB488R

GB495R

GB496R

GB497R

GB690R



GB256R

ECCOS SET

for Acculan 4 dermatome

 106 mm


consisting parts:

GB487R

GB488R

GB498R



















JC254R



AESCULAP® MECHANICAL PROCESSING

ECCOS

Acculan 4

<p>GB482R</p> <p>ECCOS HOLDER</p> <p>for Acculan 4 small drill GA344</p>  106 mm	<p>GB715R</p> <p>ECCOS HOLDER</p> <p>for 1 attachment GB870R, GB871R, GB872R, GB878R, GB881R, GB884R, GB892R, GB896R, GB897R</p>  106 mm	<p>GB716R</p> <p>ECCOS HOLDER</p> <p>for 3 attachments GB870R, GB871R, GB872R, GB878R, GB881R, GB884R, GB892R, GB896R, GB897R</p>  106 mm
 	 	 
<p>GB495R</p> <p>ECCOS HOLDER</p> <p>for GA330, GA331, GA334</p>  106 mm	<p>GB496R</p> <p>ECCOS HOLDER</p> <p>for GB620R, GB621R, GB623R, GB639R, GB641R, GB645R, GB660R, GB663R, GB664R, GB665R, GB657R, GB667R, GB668R, GB669R, GB670R</p>  106 mm	<p>GB497R</p> <p>ECCOS HOLDER</p> <p>for GB620R, GB621R, GB623R, GB639R, GB641R, GB645R, GB660R, GB663R, GB664R, GB665R, GB657R, GB667R, GB668R, GB669R, GB670R</p>  106 mm
 	 	 

GB487R**ECCOS HOLDER**

for GA675

 56 mm
**GB488R****ECCOS HOLDER**


for GA678

 76 mm
**GB489R****ECCOS HOLDER**


for GA679

 106 mm
**GB498R****ECCOS HOLDER**

for GA336, GA337, GA340, GA341

 106 mm
**GB686R****ECCOS HOLDER**















for GB436R

 56 mm


AESULAP[®] MECHANICAL PROCESSING

ECCOS HOLDER

Acculan 4

GB485R ECCOS HOLDER for GA875  56 mm	GB486R ECCOS HOLDER for GA878  56 mm	GB688R ECCOS HOLDER for BA720R  106 mm
 	 	 
GB481R ECCOS HOLDER for saw blades and other accessories  56 mm	GB690R ECCOS HOLDER for keys GA672207, GA344244	
 	 	


ECCOS HOLDER

Universal Tools


GB459R	JG303R	JG340 / JG342 / JG344
ECCOS® MOUNTING KIT	REPLACEMENT SCREWS-SET	SILICONE CLAMP
5 mm / 8.5 mm / 12 mm		
		

AESCULAP® POWER SYSTEMS


TOOL ATTACHMENTS OVERVIEW




AO LARGE




AO SMALL




DENTAL




DIN




HARRIS




HUDSON




HEXAGON



SYNTHES
DHS/DCS
















TRINKLE



ZIMMER

SIGNS & SYMBOLS


	<p>LATEX-FREE</p> <p>All products listed in this catalogue are latex-free</p>
	<p>ARTICLE CODE ECCOS FIXATION</p>
	<p>FORWARD / REVERSE DRIVE</p>
	<p>FORWARD DRIVE</p>
	<p>OSCILLATION FREQUENCY</p>
	<p>OSCILLATON ANGLE</p>
	<p>STROKE</p>
	<p>CANNULATION</p>
	<p>CABLE LENGTH</p>
	<p>WEIGHT</p>
	<p>PACKAGING UNIT</p>
	<p>SINGLE USE</p>
	<p>REQUIRED HEIGHT OF TRAY</p>
	<p>STERILE</p>

AESCULAP® POWER SYSTEMS

BROCHURES

Main Catalogues

084002 Power Systems



017599 Burrs & Blades




ELAN 4


071602 ELAN 4



071511 Mechanical Cleaning



012002 Manual Cleaning




Acculan 4


001802 Acculan 4




001902 Acculan 4 Dermotome




002002 Acculan 4 Electro





002102 Acculan 4 Sternum Saws




INDEX

Art. No.	
BA720R	40
GA031R	50
GA062R	50
GA258SU	31
GA259SU	31
GA261SU	31
GA262SU	31
GA263SU	31
GA264SU	31
GA265SU	31
GA266SU	31
GA267SU	31
GA268SU	31
GA269SU	31
GA320	46
GA330	34
GA331	38
GA334	39
GA336	38
GA337	38
GA340	40
GA341	40
GA344	42
GA346	46
GA395SU	31
GA675	46
GA676	46
GA678	46, 49
GA679	46
GA702R	21
GA703R	21
GA705R	22
GA706R	22
GA707R	22
GA708	22
GA710R	21
GA712R	21
GA713R	21
GA715R	21
GA722	22
GA724	29
GA731	29
GA732	29
GA733	29
GA736	29
GA737	29
GA749	23
GA760	27
GA761	24
GA762	24
GA763	24
GA764	24
GA765	25

Art. No.	
GA766	25
GA767	25
GA768	25
GA769	25
GA800	7
GA804	8
GA805	8
GA806	7
GA808	7
GA810	7
GA822	8
GA824	8
GA831	15
GA832	15
GA833	15
GA836	15
GA837	15
GA844	16
GA849	9
GA860	13
GA861	10
GA862	10
GA863	10
GA864	10
GA865	11
GA866	11
GA867	11
GA868	11
GA869	11
GA875	49
GA876	49
GA877	49
GA878	49
GB053R	58
GB056R	58
GB057R	58
GB062R	58
GB063R	58
GB064R	59
GB065R	58
GB066R	58
GB067R	58
GB068R	58
GB072R	60
GB073R	55
GB076R	55
GB077R	55
GB082R	55
GB083R	55
GB084R	57
GB085R	55
GB086R	55
GB087R	55

Art. No.	
GB088R	55
GB176R	50
GB184R	50
GB243R	61
GB256R	61
GB257R	61
GB260R	54
GB261R	54
GB262R	54
GB436R	39
GB459R	65
GB481R	64
GB482R	56, 62
GB485R	64
GB486R	64
GB487R	63
GB488R	63
GB489R	63
GB495R	62
GB496R	62
GB497R	62
GB498R	63
GB572R	60
GB600	52
GB600840	52
GB600850	52
GB600860	52
GB600870	52
GB600880	52
GB620R	34
GB621R	34
GB623R	34
GB628R	50
GB639R	35
GB641R	37
GB645R	35
GB657R	36
GB660R	37
GB663R	35
GB664R	35
GB665R	35
GB667R	36
GB668R	36
GB669R	36
GB670R	36
GB679R	57, 59
GB686R	63
GB688R	64
GB690R	64
GB691R	59
GB692R	57
GB697R	59
GB698R	57

Art. No.	
GB699R	57, 59
GB715R	56, 62
GB716R	56, 62
GB718R	60
GB719R	60
GB720R	60
GB723R	57, 59
GB727R	60
GB729R	60
GB730R	60
GB796SU	31
GB797SU	31
GB861P	30
GB862P	30
GB863P	30
GB866R	30
GB870R	16, 42
GB871R	16, 42
GB872R	16, 42
GB878R	18, 44
GB881R	17, 43
GB882R	17, 43
GB884R	17, 43
GB886R	17, 43
GB887R	17, 43
GB891R	18, 44
GB892R	18, 44
GB894R	18, 44
GB896R	19, 45
GB897R	19, 45
GB920R	13, 27
GB921R	13, 27
GB922R	13, 27
GB925R	13, 27
GB926R	13, 27
GB927R	13, 27
GB941R	9, 23
GB942R	9, 23
GB943R	9, 23
GB945R	9, 23
GB947R	9, 23
JG303R	65
JG340	65
JG342	65
JG344	65
TE730	30, 50
TE734	30, 50
TE735	30, 50
TE780	30, 50
XF457R	50

NOTES

AESCULAP® – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 7461 95-0 | Fax +49 7461 95-2600 | www.aesculap.com

The main product trademark "Aesculap" is a registered trademark of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.